

# **EuroBioSciences**

**Data Sheet** 

# anti-human CD34 FITC-conjugated

Cat-No.: H12166F 1 ml

#### Clone:

4H11[APG]

#### Specificity:

The antibody reacts with Class III epitope on CD34 antigen. The antibody 4H11[APG] completely blocks binding of Class II antibody QBEnd10 and Class III antibodies BIRMA K3 and 8G12 on KG1a cell line.

## Isotype subclass:

mouse IgG1

#### Form:

The purified antibody is conjugated with Fluoresceinisothiocyanate (FITC) under optimum conditions. The reagent is free of unconjugated FITC and adjusted for direct use. No reconstitution is necessary.

### Physical state:

Liquid

# **Buffer/Additives/Preservative:**

PBS containing BSA and 15 mM sodium azide, pH 7.4

# **Expiration date:**

The reagent is stable until the expiry date stated on the vial label

#### Storage conditions:

Store at 4 °C. Do not freeze. Avoid prolonged exposure to the light.

# Application:

Flow Cytometry

### References:

Leukocyte Typing VI. Kishimoto T. et al. (Eds.), Garland Publishing Inc. (1997)

#### Warning:

Sodium azide is harmful if swallowed (R22). Keep out of reach of children (S2). Keep away from food, drink and animal feeding stuff (S13). Wear suitable protective clothing (S36). If swallowed, seek medical advice immediately and show this container or label (S46). Contact with acids liberates very toxic gas (R32). Azide compounds should be flushed with large volumes of water during disposal to avoid deposits in lead or copper plumbing where explosive conditions can develop.

This material is offered for <u>research only</u>. Not for use in human. For in vitro use only. EuroBioSciences will not be held responsible for patent infringement or other violations that may occur with the use of our products.

Phone: +49 (0) 4491-9387804, Fax: +49 (0) 4491-9387805

E-Mail: info@eurobiosciences.com