

EuroBioSciences

Data Sheet

anti-rat CD4

Cat-No.: R32126 0.1 mg

Clone: W3/25

Specificity: This anti-rat T helper cell monoclonal antibody recognizes a derteminant on the majority of thymocytes (90-95%), a subset of peripheral T cells and peritoneal macrophages. (1,2,3,10.) The antigen recognized by this antibody is a surface glycoprotein of Mr 48 000 - 52 000 and is the homologue of the human CD4 and the mouse L3/T4 antigen. This antibody labels the rat T helper subset, which mediates the helper acivity for B and T cells, grafts vs. host (GVH) reactivity and produces IL-2 in the mices lymphocyte reaction (MLR). (2,4,6) Addition of the antibody to the MLR, inhibited proliferation and blocks the production of IL-2. (4,5,6,9) T cells which mediate cytotoxicity and suppressor functions are not labelled. (Thus, cells labelled by this antibody are not labelled by MRC OX-8.). This clone is invaluable for seperating T cells subsets for funcional studies. It has been studying the role of T lymphocytes in graft rejectons (7) and in studying the subsets of T cells in the rat which mediate graft vs. host disease. (8) This antibody is one of three antibodies which labels T lymphocyte populations in the rat. These clones include W3/13, which labels all T cells, as well as MRC OX-8 and W3/25 which label non-overlapping T cell subpopulations. These monoclonal antibodies used in concert are being employed extensively to investigate cellular aspects of the immune response in rats and prove to be useful as markers of functionally distinct subpopulations of lymphocytes.

Isotype subclass: Mouse IgG 1

Form: purified

Physical state: Liquid

Buffer/Additives/Preservative: PBS containing 0.09 % sodium azide (pH 7.4).

Expiration date: The reagent is stable until the expiry date stated on the vial label.

Storage conditions: For long term storage, aliquot and freeze unused portion at -20°C in volumes appropriate for single usage. Avoid freeze/thaw cycles.

Application:

Flow Cytometry

Immunhistochemistry - The clone W3/25 is reported to work with frozen and paraffin sections (11).

- 1. Williams, A.F., Galfre, G. and C. Milstein. (1977) Cell, 12, 6330673.
- 2. Bredeau, R.J., Carter, P.B., McMaster, W.R., Mason, D.W. and A.F. Williams. (1980) Eur. J of Immunol., 10, 609-615.
- 3. Barclay, A.N. (1981) Immunology, 42, 593-600.
- 4. Cantrell, D.A., Robins, R.A. and R.W. Baldwin. (1982) Cell Immunol. 70, 367-372.
- 5. Mason, D.W., Pugh, C.W. and M. Webb. (1981) Immunology, 44, 75-87.
- 7. Dallman, M.J., Mason, D.W. and M. Webb. (1982) Eur. J. Immunol. 12, 511-518
- 8. Mason, D.W. (1981) Transplantation, 32, 222-226.
- 9. White, R.A.H., Mason, D.W., Williams, A.F., Galfre, G. and C. Milstein. (1978) J. Exp. Med. 148, 644-673 10. Jefferies, W.A., Green, J.R. and A.F. Williams. (1985) J. Exp. Med. 162, 117-127. 11. Whiteland, J.L. et al (1995). J.Histochem. Cyochem. 43: 313-320

Sodium azide is harmful if swallowed (R22). Keep out of reach of children (S2). Keep away from food, drink, and animal feeding stuff (S13). Wear suitable protective clothing (S36). If swallowed, seek medical advice immediately and show this container or label (S46). Contact with acids liberates very toxic gas (R32). Azide compounds should be flushed with large volumes of water during disposal to avoid deposits in lead or copper plumbing where explosive conditions can develop.

This material is offered for research only. Not for use in human. For in vitro use only. EuroBioSciences will not be held responsible for patent infringement or other violations that may occur with the use of our products.

Phone: +49 (0) 4491-9387804, Fax: +49 (0) 4491-9387805

E-Mail: info@eurobiosciences.com